

E5

Class N autoclave for dental technicians



E5 is an autoclave from the Euronda Pro System sterilization line that offers an outstanding quality-price ratio while still maintaining Euronda's exceptional safety standards.

It is a completely automatic Class N electronic autoclave that is controlled by a microprocessor. It features a thermodynamic working method with pressure and temperature control systems. There are two sterilization cycles (121°C and 134°C), allowing users to select the one that best suits the needs of their dental practice. Their progress can be tracked thanks to the light and sound signals from the control panel. Drying time: 15 minutes. Door opening/closing mechanism with dual safety microswitch. Anti-reflow separated tanks, analogical pressure gauge, synoptic cycle advancement.



Classification:

Medical Device Class IIb
Moulded chamber
in stainless steel and
certified to PED 97/23/EC
(Pressure Equipment
Directive)



Chamber capacity:

18 L



Weight:

35 Kg (machine empty
with basket and trays)



External dimensions (wxhxd):

450 x 385 x 500 mm



Power consumption:

1500 W 6.5 A