

Euronda's range of stools includes three models designed to satisfy various needs. CDS301, Coral and Onyx are perfect for those who spend a long time sitting down thanks to their curved backrest, comfortable seat and sturdy structure.

They're also very stylish and come in a vast range of colours in both stitched and thermoformed versions.

You're spoilt for choice.

Heat fused seats



Seamed seats



The colours pictured here are purely representative.

Euronda®

Over 40 years of protection in dentistry

Euronda, synonymous with outstanding protection in the dental field for over 40 years. With an Italian origin and an international vocation, it has 4 branches operative and sells products in over 100 countries worldwide. Euronda is one of the benchmark companies for infection control in the dental sector. Production includes two specific product lines, plus Euronda Alle® sterile disposable items. Euronda Monoart®, a wide range of disposable products, accessories and personal protective equipment that teams maximum product quality with vibrant colours. Euronda Pro System, the Euronda Total Protection System. A complete line of products to optimise the whole sterilisation process and make it safer.

All Euronda products guarantee maximum quality and protection standards for dentists, professionals and patients.

Euronda Spa

Via Chizzalunga 1
36066 Sandrigo
Vicenza, Italy
t (+39) 0444 656111
f (+39) 0444 656199
m info@euronda.com
www.euronda.com

España

Euronda Spa
Identyd, Av. Cabezo
de Torres 4, 30007 Murcia,
España
t 900102034
m info@euronda.es
www.euronda.es

Russia

OOO Finstar
5/2a Varshavskaya str.
196128 St. Petersburg,
Russia
t +7 (812) 635 88 94
m info@finstartech.ru
www.euronda.com.ru

Deutschland

Euronda Deutschland GMBH
Am Landwehrbach 5
48341 Altenberge,
Deutschland
t (+49) 2505 9389 0
f (+49) 2505 9389 29
m info@euronda.de
www.euronda.de

France

Euronda France Eurl
ZAC les Vallées, Avenue de
Bruxelles, 60110 Amblainville,
France
t +33 (0) 3 44 06 69 70
f +33 (0) 3 44 02 03 89
m info@euronda.fr
www.euronda.fr

Euronda®

Operator's stool

Euronda®
Studio

CDS301

Robust, stable and safe to sit on for lengthy periods

Chromium-plated steel structure operator's stool with five spoke base and fitted with five self-pivoting wheels. Operating mechanism incorporated in the lower part of the seat that can be 360° operated by both hands without looking for any control. Solid enveloping anatomical backrest, adjustable in three different positions. It allows a relaxed working position that minimizes the load on the back, pelvis and legs. Smooth rounded surface to ensure maximum hygiene and cleanliness.



Colours:
25 colours available



Models:
Seamed seats or heat fused seats



Dimensions:
Round seat with a diameter of:
40 x 10 cm (ØxH)
Adjustable height:
51 – 69 cm



Weight:
15 kg



Coral

A practical, reliable and durable stool

Operator's stool in full compliance with ISO 7493. Structure in recyclable extruded polyamide, fireproof, scratch resistant and knock-proof, fitted with five self-locking wheels. Operating mechanism covered by a special casing and incorporated in the lower part of the seat that can be 360° operated by both hands without looking for any control. Solid enveloping anatomical backrest that can be adjusted to three positions. The seat is made of non-deformable polyurethane and covered with sky.



Colours:
25 colours available



Models:
Seamed seats or heat fused seats



Dimensions:
Round seat with a diameter of:
39 x 9 cm (ØxH)
Adjustable height:
52 – 70 cm



Weight:
13 kg



Onyx

Comfort and outstanding quality for all

Operator's stool made of an aluminium structure, fitted with 5 wheels. It's equipped with an hand lever situated under the seat for an easier regulation in height, that can be 360° operated. Round seat. Enveloping anatomical backrest, adjustable in three positions (right, centre and left) and in depth.



Colours:
25 colours available



Models:
Seamed seats or heat fused seats



Dimensions:
Round seat with a diameter of:
39 x 9 cm (ØxH)
Base diameter: 58 cm
Adjustable height:
50 – 68 cm



Weight:
11 kg

